0 – Equipment and tools needed for mounting the ceiling column

Note:

Please do refer to the “Installation & Service Manual” of the generator system prior to the mechanical and electrical installation of the ceiling mount.

All equipment and tools described in “Installation & Service Manual” Section 4 § 1.

Set of flat and pipe wrenches.

2 strings from minimum 2 meters or 80” length (or 3 strings in case of lightning system to be installed).

1 – Mechanical specifications

The ceiling plate, which supports the radiology unit, must be securely fixed to the ceiling. Choose the attachment parts that are suitable for the ceiling characteristics and that can resist a breakout force of:

- 240 Kg (529 lb.) per anchor point if there are 4 anchor points (positions shown in **Fig.1**)
- 100 Kg (220 lb.) per anchor point if there are 8 anchor points.

The plugs for heavy attachments, which are supplied with the appliance, must be used if the ceiling is made of concrete and has a minimal resistance of 23 Mpa. The plugs supplied meet the requirements specified above.

The secure and correct attachment to the ceiling must respect urban planning department regulations and must be ensured by the specialist carrying out the work. Only attachment parts that are approved by urban planning departments (plugs for heavy attachments) may be used to fit a radiology unit to the ceiling.

2 – Electrical specifications

Refer to section 4 chapter 3 of the “Installation & Service Manual”. The line cable, the control timer cable, the warning lamp cable and the RVG link cable (if used), must reach the ceiling plate.

3 – Installation

3.1 - Fitting the plate (Fig. 1)

Place the ceiling plate in the required position. Mark the 4 or 8 attachment points to the ceiling. Drill the holes into the ceiling and install the fitting system that is suitable for the nature of the ceiling. Place two strings (S) attached to the ceiling plate as in Fig. 2 and secure the plate to the ceiling before locking it

Note:

In case the column supports a lightning system, place one extra string on the plate in order to pull up the light supply cable.

3.2 - Mounting the column (Fig. 2, 3)

Introduce the strings inside the column and insert it on the neck of the plate. Lock it in place by using the security screw (a) and the nut (b).

Position the 8 pressure screws (c) around the neck of the ceiling plate until they are in contact. Adjust them until the column is vertical.

IMPORTANT:

When installation is finished, rotation of arm around the column shows a dead zone (A) that has to be positioned initially on the opposite to the working area (B). The 4 holes for screws at the bottom of the column are defining this dead zone, as shown on figure 4. Position properly the column by choosing the correct orientation before locking it to the plate (screws opposite to working area).

3.3 - Fitting the power board (Fig. 5, 6, 7)

Take the board support plate and identify the face towards the power board (sticker on it). Install on it as indicated on figure 5 the spacers and pins: 4 plastic pins (A) and 2 metal ground pins (B).

Fix the board support plate to the ceiling plate using 3 screws and washers (C) as in fig.6. Then place the insulation sheet.

Take the power board out of its box. In the case of Kodak 2200, do not forget to install the adaptor board on it. Then fit the board using the nuts supplied on the spacers mounted on the support plate (see fig.7).

3.4 - Fitting the cover (Fig. 8)

Slide the cover over the length of the column and hold it half-way with the O-ring.

3.5 - Fitting the extension and the radiology unit (Fig. 9 to Fig 13 and Fig.15)

Locate the string on the line voltage input side of the power board, attach the ON/OFF switch cable end to the string (fig.9) and while pulling it approach the end of the extension arm to the column as shown in fig.10. Pass the cable between the ceiling and the plate, and through the dedicated hole in the plate (fig.13).

Introduce the other string through the bearing assembly and leave it hanging.

Introduce the bearing assembly inside the column and secure it with the 4 countersunk screws (fig.11), the screw covers, and the plastic plug on the big hole of the column.

Install the scissor arm (closed with its strap) as described in figure 15 and in section 5, chapter 2.3 of the "Installation & Service Manuel".

Pass the arm cable through the extension arm. Attach the string to its end (fig.12) and pull the cable inside the column up to the plate (fig.13).

Install the generator as described in section 5 chapter 2.4 (fig.15).

Install the brake on the arm as described in section 5 chapter 2.5.

Put the plastic covers in place.

3.6 - Electrical connections on the plate (Fig. 14)

Note

For electrical connections, see also Section 5 chapter 3.2.2 (for Kodak 2100) or chapter 4.2.2 (for Kodak 2200): "Wiring"

Connect the On/Off switch cable to the MAIN PWR connector (4).

Connect the generator cable to the PWR ARM and CTRL ARM connector (5).

Connect the cable from the power supply to the MAIN PWR connector (3).

Connect the control timer cable to the CTRL PANEL connector J10 in case of Kodak 2100 system, or on the adaptor board in the case of Kodak 2200 system (2).

Connect, if needed, the RVG synchronization cable (1).

Connect, if used, the "ready" light cable to J9 connector (7).

IMPORTANT

In the case of Kodak 2200, do not forget to close pins 2 & 3 of J7 connector with a jumper. If not done, then an error message ("POWER ERROR") will appear and exposing will not be possible.

Attach all cables to the plate using nylon cable ties.

Check the earth recovery connection.

3.7 – End of installation

Place the cover against the ceiling and hold it in place with the O-ring.

Install and connect the control timer as described in section 5:

- For Kodak 2100, chapter 2.7.1. Refer to chapter 2.7.3 in case of double exposure switches (USA) for the wiring and jumper configurations.
- For Kodak 2200, chapter 2.7.2. Refer to chapter 2.7.4 in case of double exposure switches (USA) for the wiring and jumper configurations.

4 – Adaptation kit for lightning system

There is possibility to use the ceiling column as a support for a lightning system.

With the standard configuration, there is a standard grip (fig.17 & 18) that is supplied as accessory (Cat. 5155825) that allows adapting easily most types of existing systems (fig.19).

With the Faro ceiling configuration, the grip (fig.19 & 20) is supplied in the global kit.

WARNING

Equipments fixed on the column are limited in terms of dimensions and weight (Total 15 kg / 33 lb at 1,5 meter / 59" span). On the same way, their electrical installation must comply with all current installation standards.

This part of installation is done under the responsibility of the fitter. Eastman Kodak Company will not be held responsible for any modification to the original equipment by a non-Kodak certified technician.

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